

AASRAA FLYTE GO PLUS

LIGHTWEIGHT TRAVEL WHEELCHAIR

AR117S

The AASRAA Flyte Go Plus AR117S is a **lightweight travel** wheelchair designed for convenience and durability. Weighing only **9 kg**, it features a **strong powder-coated mild steel frame**, compact **front-folding design**, and a **pull-out travel handle** for easy transport. **Honeycomb shock-absorbing tyres**, **breathable antimicrobial mesh cushioning**, **swing-up armrests**, and **dual attendant brakes** ensure comfort, safety, and smooth mobility across surfaces.

★ **MRP: INR 16999**



THE MOST LIGHTWEIGHT WHEELCHAIR

-  **100 KG**
MAX LOAD CAPACITY
-  **BREATHABLE MESH**
SEAT FABRIC TYPE
-  **460MM (18")**
SEAT WIDTH
-  **MAG WHEELS**
REAR WHEEL TYPE
-  **9 KG**
NET WEIGHT
-  **10.8 KG**
GROSS WEIGHT
-  **83 X 60 X 93 CM**
OPEN DIMENSIONS
-  **60 X 36 X 78 CM**
BOX DIMENSIONS

TRUCK SIZE	APPROX QTY
8'	12-16
14'	40-50
17'	70-80
22'	100-110
27'	130-140
34'	250-270
40'	290-310

SCAN FOR
PRODUCT
VIDEO &
DETAILED
SPECS



LOADING CAN VARY DEPENDING ON VEHICLE INNER DIMENSIONS



ULTRA COMPACT & EASY FOR TRAVEL

HIGHLIGHTS



SOLID HONEYCOMB WHEELS FOR BETTER SHOCK ABSORBANCY



ULTRA LIGHTWEIGHT AND QUICK FOLDING



TRAVEL FRIENDLY DESIGN



360° FRONT WHEEL WITH HONEYCOMB AND SHOCK ABSORBERS

FEATURES



**SAFETY FOLDING
FRAME LOCK**



**360° FRONT WHEEL
WITH PLASTIC FORK**



FLIP UP FOOT PLATE



**SOFT, FOAM PADDED
ARMPAD**



**PLASTIC TYRES WITH
100% VIRGIN PPCP
MATERIAL**



**PUNCTURE PROOF
SOLID REAR WHEELS**



**180° HAND BRAKE
WITH
PARKING LOCK**



**ADJUSTABLE SEAT
BELT**



**COMFORTABLE
SEATING AND
STABILITY**




**LIGHT WEIGHT &
EASY TO CARRY**




**TRAVEL FRIENDLY
DESIGN**



MESH CUSHION SEAT



**EASY VECLRO
FOLDING STRAP**



**270° SWING-UP
ARMRESTS**



**COMPACT &
FOLDABLE**



**BACK POCKET TO
KEEP FILES/ FOLDERS**

SPECIFICATIONS

DIMENSIONS	MM	INCH
OPEN LENGTH	830	33
OPEN LENGTH (WITHOUT FOOTREST)	670	26
OPEN WIDTH	600	24
OPEN HEIGHT	925	36
FOLDED WIDTH	355	14
FOLDED LENGTH	600	24
FOLDED HEIGHT	775	31
SEAT WIDTH	460	18
SEAT DEPTH	370	15
BACKREST HEIGHT (FROM SEAT)	330	13
BACKREST THICKNESS	13	1
SEAT THICKNESS	30	1
EXTRA CUSHION THICKNESS	N/A	N/A
ARMREST HEIGHT FROM SEAT	230	9
SEAT HEIGHT FROM FLOOR	515	20
SEAT HEIGHT FROM FOOTREST	315	12
GROUND CLEARANCE	115	5
FOOTREST HEIGHT FROM FLOOR	195	8
FRONT WHEEL DIAMETER	148	6
REAR WHEEL DIAMETER	295	12
DISTANCE BETWEEN SEAT & FOOTREST	315	12
BOX LENGTH	600	24
BOX WIDTH	360	14
BOX HEIGHT	780	31

DIMENSION	MEASUREMENT
WEIGHT OF WHEELCHAIR (NET WEIGHT)	9 KG
WEIGHT WITH PACKAGING (GROSS WEIGHT)	10.8 KG
TURNING RADIUS	1.5 M
LOAD CAPACITY	100 KG

TECHNICAL SPECIFICATIONS	
PROPULSION SYSTEM	MANUAL-ATTENDANT
FRAME MATERIAL	Powder Coated Mild Steel & Aluminium Alloy
FRAME COLOR	Silver Glossy
CROSSBAR TYPE	CROSS-BRACE SYSTEM
CROSSBAR MATERIAL	Powder Coated Mild Steel & Aluminium Alloy
CROSSBAR COLOR	BLACK MATT
SEAT TYPE	Screw Hidden + Mesh Cushion + Hard Bottom Seat
UPHOLSTREY	Breathable Mesh Cushion, Antimicrobial and Fire-Retardant
UPHOLSTREY COLOR	Black and Orange
FRONT WHEEL MATERIAL	6 Inch MAG Design PP Castor With Honeycomb PVC Tyres & Shock Absorbers
REAR WHEEL MATERIAL	12 Inch MAG Design PP Castor With Honeycomb PVC Tyres
REAR TYRE MATERIAL	PVC Solid Moulded Tyres
WHEELCHAIR TYPE	Front Compact Fold
FOOTREST TYPE	Swing Open/Close
FOOTREST MATERIAL	Aluminium + MS Alloy
LEGREST TYPE	Front Lift
BACKPOCKET	YES
ARMREST TYPE	270° SWINGABLE
ARMREST MATERIAL	FOAM PADDED
HEADREST	NO
ANTI TIPPER WHEELS	NO
PROPELLING RING	N/A
BRAKING SYSTEM	Dual Locking & Push-To-Brake Attendant Brakes - 180 Degree Adjustable
ATTENDANT BRAKES	Yes
PUSH HANDLES	YES, WITH FOAM GRIPS
PROVISION OF COMMODE	NO
CALF STRAP MATERIAL	N/A
SIDE GAURD MATERIAL	N/A
ACCESSORIES	Pull-Out Handle For Travel Purposes, Tool Kit

SUITABLE FOR USERS



SENIORS



CHILDREN



ADULTS



WOMEN



~~ATHLETES~~



~~PwD~~



~~PATIENTS~~



PREGNANT



SHORT
TERM
USERS



~~CEREBRAL
PALS~~

SUITABLE FOR ACTIVITY



FLIGHT/TRAIN



CAR



BATHROOM



OFFICE



MARKET




OUTDOOR


SUITABLE FOR SURFACE



TILE




WOOD



PAVERS




ASPHALT



CARPET



~~SAND~~



GRAVEL



~~GRASS~~



आपका भरोसा, *aasraa* का वादा

your trust, aasraa's assurance

Leading brand in India for wheelchairs (manual and electric) and other mobility and orthopaedic aids



Manufactured & Marketed By:

PEAAR HEALTHTECH LLP

Plot No 12, Sector 4
Faridabad 121004, Haryana (INDIA)

(+91) 828 747 6144 | (+91) 875 005 0480

contact@peaar.co.in



www.aasraa.co.in