

# AASRAA

## PORTACOMMODE

COMMODE TRANSFER CHAIR FOR WC

### ART05S

The **AASRAA Portacommode** makes bathroom access simple, safe, and dignified. Its height-adjustable design slides easily over any Western commode, while the padded U-cut seat ensures comfort and hygiene. With detachable armrests and backrest, it allows smooth transfers and 4-in-1 versatile use. Compact, rust-resistant, and easy to maneuver in narrow spaces—perfect for daily home care.

★ **MRP: INR 11,999.00**



THE MULTIFUNCTION BATHROOM CHAIR

-  **120 KG**  
MAX LOAD CAPACITY
-  **REXINE**  
SEAT FABRIC TYPE
-  **460MM (18')**  
SEAT WIDTH
-  **PVC WHEELS**  
WHEEL TYPE
-  **13 KG**  
NET WEIGHT
-  **15.5 KG**  
GROSS WEIGHT
-  **96X 54 X 91 CM**  
OPEN DIMENSIONS
-  **75 X 58 X 73 CM**  
BOX DIMENSIONS

| TRUCK SIZE | APPROX QTY |
|------------|------------|
| 8'         | 8-10       |
| 14'        | 30-40      |
| 17'        | 60-70      |
| 22'        | 80-90      |
| 27'        | 100-120    |
| 34'        | 200-220    |
| 40'        | 250-260    |

LOADING CAN VARY DEPENDING ON VEHICLE INNER DIMENSIONS

SCAN FOR  
PRODUCT  
VIDEO &  
DETAILED  
SPECS




FULLY DETACHABLE AND EASY TO STORE




## HIGHLIGHTS



SLIDE OVER **ALL** WC - DESIGNER, WITH TANK, WALL MOUNTED, FLOOR MOUNTED, ETC



BACK SLIDE OUT BUCKET IN CASE OF BUCKET USAGE



FULLY DETACHABLE ARMRESTS FOR EASY PATIENT TRANSFER



4 DIFFERENT WAYS TO USE - AS PER YOUR COMFORT AND WASHROOM DESIGN

## FEATURES



**5 LEVEL HEIGHT  
ADJUSTMENT**



**4 INCH REAR WHEELS  
WITH FOOT BRAKE**



**FLIP UP FOOTPLATE**



**SOFT SOLID PU  
WATERPROOF  
ARMPAD**



**4 INCH PVC FRONT  
WHEELS**



**REMOVABLE U-CUT  
CUSHION**



**HOLLOW U-CUT  
WITHOUT FRONT BAR**



**STEEL FOOT PLATE**



**EASY DETACHABLE  
BACKREST**



**PLASTIC FLIP SWITCH  
FOR ARMREST  
REMOVAL**




**BACK REMOVABLE  
BUCKET WITH LID**



**WIDTH 21" - PASSES  
THROUGH NARROW  
DOORS**



**40" PADDED FOAM  
SEAT**



**WATERPROOF  
REXINE UPHOLSTREY**



**REMOVABLE  
ARMRESTS**



**COMPACT &  
FOLDABLE**

## SPECIFICATIONS

| DIMENSIONS                                    | MM         | INCH         |
|---|------------|--------------|
| OPEN LENGTH                                   | 955        | 37.6         |
| OPEN WIDTH                                    | 540        | 21           |
| FOLDED WIDTH                                  | -          | -            |
| HEIGHT  | 905        | 35.6         |
| SEAT WIDTH                                    | 460        | 18           |
| SEAT DEPTH                                    | 445        | 17.5         |
| BACKREST HEIGHT (FROM SEAT)                   | 360        | 14           |
| BACKREST THICKNESS                            | 7          | 0.3          |
| SEAT THICKNESS                                | 48         | 2            |
| EXTRA CUSHION THICKNESS                       | -          | -            |
| ARMREST HEIGHT FROM SEAT                      | 220        | 9            |
| SEAT HEIGHT FROM FLOOR                        | 480        | 19           |
| SEAT HEIGHT FROM FOOTREST (LOWEST TO HIGHEST) | 510 TO 580 | 20 TO 23     |
| WC CLEARANCE FROM GROUND (LOWEST TO HIGEST)   | 420 to 490 | 16.5 to 19.3 |
| SEAT HEIGHT FROM FOOTREST                     | 410        | 16           |
| ALL WHEEL DIAMETER                            | 102        | 4            |
| REAR WHEEL DIAMETER                           | 580        | 22.8         |
| DISTANCE BETWEEN SEAT & FOOTREST              | 410        | 16           |
| BOX LENGTH                                    | 750        | 30           |
| BOX WIDTH                                     | 575        | 23           |
| BOX HEIGHT                                    | 725        | 29           |

| DIMENSION                            | MEASUREMENT |
|--------------------------------------|-------------|
| WEIGHT OF WHEELCHAIR (NET WEIGHT)    | 13 KG       |
| WEIGHT WITH PACKAGING (GROSS WEIGHT) | 15.8 KG     |
| TURNING RADIUS                       | 1.00 M      |
| LOAD CAPACITY                        | 125 KG      |

| TECHNICAL SPECIFICATIONS |  |
|--------------------------|--|
| PROPULSION SYSTEM        | MANUAL-ATTENDANT   |
| FRAME MATERIAL           | Powder Coated Mild Steel (Carbon Steel) (MS ERW 22.2x1.2 MM - 18 GAUGE)<br>Powder Coated Mild Steel (Carbon Steel) (MS ERW 25.4x1.2 MM - 18 GAUGE)<br>Powder Coated Mild Steel (Carbon Steel) (MS ERW 19.05x1.2 MM - 18 GAUGE) |
| FRAME COLOR              | SILVER GLOSSY  |
| SEAT TYPE                | U-Cut Commode Seat   |
| UPHOLSTREY               | Waterproof Rexine Leatherite   |
| UPHOLSTREY COLOR         | BLACK  |
| EXTRA CUSHION            | N/A  |
| EXTRA CUSHION COLOR      | N/A  |
| FRONT WHEEL MATERIAL     | MS Castor 4 Inch with PVC Tyre Without Brake   |
| REAR WHEEL MATERIAL      | MS Castor 4 Inch with PVC Tyre With Brake  |
| REAR TYRE MATERIAL       | PVC  |
| WHEELCHAIR TYPE          | DETACHABLE   |
| FOOTREST TYPE            | Swing Open/Close   |
| FOOTREST MATERIAL        | Powder Coated Mild Steel (Carbon Steel) (MS ERW 22.2x1.2 MM - 18 GAUGE)  |
| BACKPOCKET               | NO   |
| ARMREST TYPE             | DETACHABLE   |
| ARMREST MATERIAL         | MOULDED SOFT PU  |
| HEADREST                 | NO   |
| ANTI TIPPER WHEELS       | N/A  |
| PROPELLING RING          | N/A  |
| BRAKING SYSTEM           | FOOT BREAKS IN REAR WHEEL  |
| ATTENDANT BRAKES         | NO   |
| PUSH HANDLES             | YES, WITH PVC GRIPS  |
| PROVISION OF COMMUNE     | YES  |
| CALF STRAP MATERIAL      | N/A  |
| SIDE GAURD MATERIAL      | N/A  |
| ACCESSORIES              | BUCKET WITH LID  |

## SUITABLE FOR USERS



SENIORS



CHILDREN



ADULTS



WOMEN



ATHLETES



PWD



PATIENTS



PREGNANT



SHORT  
TERM  
USERS



~~CEREBRAL  
PALSY~~

## SUITABLE FOR ACTIVITY



FLIGHT/TRAIN



CAR



BATHROOM



OFFICE




MARKET



OUTDOOR


## SUITABLE FOR SURFACE



TILE




WOOD



PAVERS



~~ASPHALT~~



CARPET



~~SAND~~



~~GRAVEL~~




~~GRASS~~



# आपका भरोसा, *aasraa* का वादा

*your trust, aasraa's assurance*

Leading brand in India for wheelchairs (manual and electric) and other mobility and orthopaedic aids



**Manufactured & Marketed By:**

**PEAAR HEALTHTECH LLP**

Plot No 12, Sector 4  
Faridabad 121004, Haryana (INDIA)

(+91) 828 747 6144 | (+91) 875 005 0480

contact@peaar.co.in



[www.aasraa.co.in](http://www.aasraa.co.in)