

AASRAA RED ROVER

BASIC ROLLATOR WALKER

ARR05R

A sturdy, easy-to-use rollator designed for everyday convenience. With **height-adjustable handles**, a **padded seat and back-support**, and **dual-action brakes (push to slow, pull to lock)**, it offers safety and support whether you're walking or resting. Equipped with a **durable metal basket** for personal items, **smooth-rolling PVC front wheels with metal forks**, and stable **dual rear wheels**, it's built to move with you – indoors or out. When not in use, it **folds easily** into a **compact size** for storage or travel.

★ **MRP: INR 13,999.00**



BEST QUALITY ROLLATOR WALKER

-  **100 KG**
MAX LOAD CAPACITY
-  **REXINE**
SEAT FABRIC TYPE
-  **800 - 900MM**
MIN TO MAX HEIGHT
-  **MILD STEEL**
FRAME MATERIAL
-  **7.9 KG**
NET WEIGHT
-  **9.2 KG**
GROSS WEIGHT
-  **67 X 52 X 80 CM**
OPEN DIMENSIONS
-  **62 X 25 X 77 CM**
BOX DIMENSIONS

TRUCK SIZE	APPROX QTY
8'	35-35
14'	100-120
17'	130-150
22'	160-180
27'	190-210
34'	370-400
40'	430-450

LOADING CAN VARY DEPENDING ON VEHICLE INNER DIMENSIONS

SCAN FOR
PRODUCT
VIDEO &
DETAILED
SPECS




LIGHTWEIGHT, COMPACT & FOLDABLE



FEATURES



**WATERPROOF SEAT
WITH 40MM SOFT
CUSHION**




**SOFT BACK SUPPORT
FOR COMFORT**



**FRONTWHEEL WITH
METAL BRACKET**



**DUAL-ACTION BRAKES
(PUSH TO SLOW, PULL
TO LOCK)**



**SOFT & COMFTABLE
PVC GRIP**



**FLIP OPEN SEAT
WITH BASKET FOR
ESSENTIALS**



**DOUBLE REAR
WHEELS WITH
BRAKES**



**LIGHTWEIGHT &
COMPACT**



**HEIGHT ADJUSTABLE
HANDLES**

SPECIFICATIONS

DIMENSIONS	MM	INCH
OPEN LENGTH (AT BOTTOM)	665	26
OPEN WIDTH (AT BOTTOM)	520	20
OPEN WIDTH (AT HANDLEBAR)	515	20
HEIGHT (LOWEST TO HIGHEST)	800-900	31-35
FOLDED LENGTH	570	22
FOLDED WIDTH	270	11
FOLDED HEIGHT	750	30
GRIP DISTANCE FROM SCREW	N/A	N/A
GRIP EXTRA FOAM THICKNESS	N/A	N/A
RECIPROCAL MAX RANGE	N/A	N/A
GAP BETWEEN EACH HEIGHT ADJUSTMENT (RANGE)	40	2
TOTAL LEVELS OF HEIGHT ADJUSTMENT NUMBER	4	
GROUND CLEARANCE	210	8
SEAT WIDTH	400	16
SEAT DEPTH	300	12
SEAT HEIGHT FROM FLOOR	500	20
BACKREST HEIGHT (FROM SEAT)	380	15
SEAT THICKNESS	40	2
HANDLE HEIGHT FROM SEAT	300-400	12-16
FRONT WHEEL DIAMETER	150	6
REAR WHEEL DIAMETER	150	6
HANDLE HEIGHT FROM GROUND (LOWEST TO HIGHEST)	800-900	31-35
BOX LENGTH	615	24
BOX WIDTH	250	10
BOX HEIGHT	765	30

TECHNICAL SPECIFICATIONS	
ASSEMBLY SYSTEM	SEMI-ASSEMBLED
FRAME MATERIAL	Powder Coated Carbon Steel (MS ERW 25.4x1.2 MM - 18 GAUGE)
FRAME COLOR	RED GLOSSY
CENTRE BAR MATERIAL	N/A
CENTRE BAR COLOR	N/A
SEAT TYPE	Soft Cushion Padded Rexine Seat
BASKET TYPE	Metal Cage
BASKET COLOR	BLACK
UPHOLSTREY	WATERPROOF REXINE, Antimicrobial and Fire-Retardant
UPHOLSTREY COLOR	BLACK
FOLDABLE / FIXED	FOLDABLE
PUSH HANDLES	YES, WITH PVC Soft Grip
EXTRA GRIP COVER	N/A
EXTRA CUSHION	WASHABLE CUSHION
ADJUSTMENT LEVEL	Height Adjustable Upto 4 Levels
FOOT TIP MATERIAL	N/A
FRONT WHEEL MATERIAL	Single 6" MAG Design PP Castor - With PVC Tyre And MS Fork
REAR WHEEL MATERIAL	Dual 6" MAG Design PP Castor - With PVC Tyre
BRAKING SYSTEM	Dual Mode - Push To Apply and Pull to Lock
RECIPROCAL	NO
PROVISION OF COMMODE	NO
ACCESSORIES	BASKET

DIMENSION	MEASUREMENT
WEIGHT OF PRODUCT (NET WEIGHT)	7.9 KG
WEIGHT WITH PACKAGING (GROSS WEIGHT)	9.2 KG
TURNING RADIUS	1 M
LOAD CAPACITY	100 KG

SUITABLE FOR USERS



SENIORS



CHILDREN



ADULTS



WOMEN



ATHLETES



PWD



PATIENTS



~~PREGNANT~~



SHORT
TERM
USERS



~~CEREBRAL
PALSY~~

SUITABLE FOR ACTIVITY



FLIGHT/TRAIN



CAR



BATHROOM



OFFICE




MARKET



OUTDOOR


SUITABLE FOR SURFACE



TILE




WOOD



PAVERS



ASPHALT



CARPET



~~SAND~~



~~GRAVEL~~




GRASS



आपका भरोसा, *aasraa* का वादा

your trust, aasraa's assurance

Leading brand in India for wheelchairs (manual and electric) and other mobility and orthopaedic aids



Manufactured & Marketed By:

PEAAR HEALTHTECH LLP

Plot No 12, Sector 4
Faridabad 121004, Haryana (INDIA)

(+91) 828 747 6144 | (+91) 875 005 0480

contact@peaar.co.in



www.aasraa.co.in