

AASRAA STURDI

DELUXE FOLDABLE MANUAL WHEELCHAIR

AR803E

Built for **strength, comfort, and everyday ease**. Its **foldable** design with **detachable legrests** makes **travel** simple, while **swing-up desk armrests** allow **easy transfers** and **comfortable seating at tables**. Strong **MAG wheels** with durable **rubber tyres** ensure smooth rides indoors or outdoors. With **padded seats, safety belt, pouch, and back pocket**, it's a reliable wheelchair for **daily living and travel**.

★ **MRP: INR 13,999.00**



THE MOST PREMIUM WHEELCHAIR



120 KG
MAX LOAD CAPACITY



NYLON CANVAS
SEAT FABRIC TYPE



450MM (18")
SEAT WIDTH



MAG WHEELS
REAR WHEEL TYPE



16.5 KG
NET WEIGHT



18.5 KG
GROSS WEIGHT



105 X 64 X 90 CM
OPEN DIMENSIONS



82 X 29 X 92 CM
BOX DIMENSIONS

TRUCK SIZE	APPROX QTY
8'	12-16
14'	40-50
17'	70-80
22'	100-110
27'	130-140
34'	250-270
40'	290-310

LOADING CAN VARY DEPENDING
ON VEHICLE INNER DIMENSIONS

SCAN FOR
PRODUCT
VIDEO &
DETAILED
SPECS



FEATURE RICH UNLIKE ANYTHING ELSE

HIGHLIGHTS



SWING UP ARMRESTS FOR PATIENT TRANSFERS



SCREW HIDDEN SEAT WITH WASHABLE CUSHION



DETACHABLE LEGREST TO MAKE IT COMPACT FOR TRAVEL



DESK ARMREST DESIGN TO MAKE SITTING AT TABLE COMFORTABLE

FEATURES



**HEIGHT ADJUSTABLE
FOOTREST WITH
PUSH BUTTON**




**360° FRONT WHEEL
WITH METAL FORK**



FLIP UP FOOTPLATE



**SOFT, PADDED
ARMPAD**



**RUBBER TYRES WITH
GRIP**



**SOLID MAG REAR
WHEELS**



PARKING BRAKES



**ADJUSTABLE SEAT
BELT**



**FLAT-TOP END-CAP
FOR TILT**



**LIGHT WEIGHT &
EASY TO CARRY**



SELF PROPELLING



BACK POCKET




**ANTI SLIP PLASTIC
FOOTREST**



**DETACHABLE CALF
GAURD WITH PIPING**



**COMPACT &
FOLDABLE**



**POUCH TO KEEP
PHONE/WALLET**

SPECIFICATIONS

DIMENSIONS	MM	INCH
OPEN LENGTH	1050	42
OPEN WIDTH	640	25
FOLDED WIDTH	270	11
HEIGHT	905	35.6
SEAT WIDTH	475	18.7
SEAT DEPTH	400	16
BACKREST HEIGHT (FROM SEAT)	380	15
BACKREST THICKNESS	7	0.3
SEAT THICKNESS	7	0.3
EXTRA CUSHION THICKNESS	5.5	0.2
ARMREST HEIGHT FROM SEAT	270	10.5
SEAT HEIGHT FROM FLOOR	480	19
SEAT HEIGHT FROM FOOTREST (LOWEST TO HIGHEST)	410 TO 505	16 TO 20
GROUND CLEARANCE	120	4.7
FOOTREST HEIGHT FROM FLOOR (LOWEST TO HIGHEST)	25 TO 120	1 TO 5
FRONT WHEEL DIAMETER	193	8
REAR WHEEL DIAMETER	580	22.8
DISTANCE BETWEEN SEAT & FOOTREST	400	16
BOX LENGTH	820	32
BOX WIDTH	290	11
BOX HEIGHT	920	36

DIMENSION	MEASUREMENT
WEIGHT OF WHEELCHAIR (NET WEIGHT)	16.5 KG
WEIGHT WITH PACKAGING (GROSS WEIGHT)	18.5 KG
TURNING RADIUS	1.65 M
LOAD CAPACITY	120 KG

TECHNICAL SPECIFICATIONS	
PROPULSION SYSTEM	MANUAL-SELF
FRAME MATERIAL	Powder Coated Mild Steel (Carbon Steel) (MS ERW 22.2x1.2 MM - 18 GAUGE)
FRAME COLOR	BLUE GLOSSY
CROSSBAR TYPE	DOUBLE SCISSOR PIPE MECHANISM
CROSSBAR MATERIAL	Powder Coated Mild Steel (Carbon Steel) (MS ERW 25.4x1.5 MM - 16 GAUGE)
CROSSBAR COLOR	BLACK GLOSSY
SEAT TYPE	SCREW HIDDEN
UPHOLSTREY	Nylon Canvas, Antimicrobial and Fire-Retardant
UPHOLSTREY COLOR	BLACK
EXTRA CUSHION	WASHABLE CUSHION
EXTRA CUSHION COLOR	BLUE
FRONT WHEEL MATERIAL	MAG Design PP Castor - 6 Spoke With PVC Tyre
REAR WHEEL MATERIAL	MAG PPCP Wheel - 6 Spoke
REAR TYRE MATERIAL	SOLID NATURAL RUBBER
WHEELCHAIR TYPE	FOLDABLE
FOOTREST TYPE	Height Adjustable, Swing Open/Close
FOOTREST MATERIAL	POLYMER PPCP
LEGREST TYPE	DETACHABLE
BACKPOCKET	YES
ARMREST TYPE	SWINGABLE (FLIP UP)
ARMREST MATERIAL	FOAM PADDED REXINE
HEADREST	NO
ANTI TIPPER WHEELS	NO
PROPELLING RING	SMOOTH POLYMER PPCP DESIGN
BRAKING SYSTEM	ANTI-SKID MANUAL PARKING BRAKE SYSTEM
ATTENDANT BRAKES	NO
PUSH HANDLES	YES, WITH PVC GRIPS
PROVISION OF COMMODE	NO
CALF STRAP MATERIAL	YES, NYLON CANVAS WITH PIPING
SIDE GAURD MATERIAL	POLYMER PPCP
ACCESSORIES	POUCH

SUITABLE FOR USERS




SENIORS



~~CHILDREN~~



ADULTS



WOMEN



~~ATHLETES~~



PWD



PATIENTS



PREGNANT



SHORT
TERM
USERS




~~CEREBRAL
PALSY~~

SUITABLE FOR ACTIVITY




FLIGHT/TRAIN



CAR



~~BATHROOM~~



OFFICE



MARKET




OUTDOOR


SUITABLE FOR SURFACE



TILE




WOOD



PAVERS



ASPHALT



CARPET



~~SAND~~



~~GRAVEL~~




GRASS



आपका भरोसा, *aasraa* का वादा

your trust, aasraa's assurance

Leading brand in India for wheelchairs (manual and electric) and other mobility and orthopaedic aids



Manufactured & Marketed By:

PEAAR HEALTHTECH LLP

Plot No 12, Sector 4
Faridabad 121004, Haryana (INDIA)

(+91) 828 747 6144 | (+91) 875 005 0480

contact@peaar.co.in



www.aasraa.co.in