

AASRAA EASYSTEP

LIGHTWEIGHT RECIPROCAL WALKERS

ARW15C / ARW25A

AASRAA Aluminium Walker is a **lightweight, foldable** walker with **fixed** and **reciprocal walking** modes. Made from strong **anodised aluminium** and available in **shiny** and **matt finish**, it offers **8-level height adjustment**, **durable PVC grips**, and **anti-slip rubber shoes** with metal washers for long life. Safe, sturdy, and ideal for daily mobility and recovery.

★ **MRP: INR 3,499.00**

AVAILABLE IN COLORS:



ARW25C
SHINY FINISH



ARW15A
MATT FINISH



ALUMINIUM WALKER CHROME FINISH
AASRAA CHROMEMATE

- 120 KG**
MAX LOAD CAPACITY
- DUAL MODES**
RECIPROCATING
- 8-LEVELS**
HEIGHT ADJUSTABLE
- PVC**
GRIP TYPE
- 2.0 KG**
NET WEIGHT
- 3.5 KG**
GROSS WEIGHT
- 53 X 45 X 74 CM**
OPEN DIMENSIONS
- 61 X 12 X 77 CM**
BOX DIMENSIONS



ALUMINIUM WALKER MAT FINISH
AASRAA EASYSTEP

TRUCK SIZE	APPROX QTY
8'	12-16
14'	40-50
17'	70-80
22'	100-110
27'	130-140
34'	250-270
40'	290-310

LOADING CAN VARY DEPENDING
ON VEHICLE INNER DIMENSIONS

SCAN FOR
PRODUCT
VIDEO &
DETAILED
SPECS



HIGHLIGHTS



SCREW WITH SAFETY CAPS




120KG LOAD CAPACITY



EASY FOLDABLE & COMPACT



LIGHTWEIGHT & EASY TO CARRY



WATERPROOF PVC GRIPS



8-LEVEL HEIGHT ADJUSTABLE



ANTI SLIP SHOES



WALKING & FIXED MODE



EASY PUSH BUTTON TO SWITCH MODE

SPECIFICATIONS

TECHNICAL SPECIFICATIONS	
ASSEMBLY SYSTEM	Fully Assembled
FRAME MATERIAL	Anodised Aluminium (18 Gauge)
FRAME COLOR	Shiny Silver Chrome / Matt Silver
FOLDABLE / FIXED	FOLDABLE
GRIP TYPE	PVC GRIP
ESTRA GRIP COVER	NO
ADJUSTMENT LEVEL	Height Adjustable Upto 8 Levels
Centre Bar Color	Metallic Grey
Foot Tip Material	Anti Slip Grip
Reciprocal	YES

DIMENSIONS	MM	INCH
OPEN LENGTH (AT BOTTOM)	530	21
OPEN WIDTH (AT BOTTOM)	450	18
TOTAL HEIGHT (LOWEST TO MAX)	740-900	29-35
GAP BETWEEN EACH HEIGHT ADJUSTMENT (RANGE)	20	1
GRIP DISTANCE FROM JOINT	40	2
RECIPROCAL RANGE	140	6
FOLDED LENGTH	530	30
FOLDED WIDTH	100	4
FOLDED HEIGHT	740	29
BOX LENGTH	610	24
BOX WIDTH	120	5
BOX HEIGHT	770	30

DIMENSION	MEASUREMENT
WEIGHT OF WALKER (NET WEIGHT)	2.1 KG
WEIGHT WITH PACKAGING (GROSS WEIGHT)	3.55 KG
TOTAL LEVELS OF HEIGHT ADJUSTMENT	8
LOAD CAPACITY	120 KG

SUITABLE FOR USERS



SENIORS



CHILDREN



ADULTS



WOMEN



ATHLETES



PWD



PATIENTS



PREGNANT



SHORT TERM USERS



CEREBRAL PALSY


SUITABLE FOR SURFACE




TILE



WOOD




PAVERS



ASPHALT



CARPET



SAND



GRAVEL



GRASS



आपका भरोसा, *aasraa* का वादा

your trust, aasraa's assurance

Leading brand in India for wheelchairs (manual and electric) and other mobility and orthopaedic aids



Manufactured & Marketed By:

PEAAR HEALTHTECH LLP

Plot No 12, Sector 4
Faridabad 121004, Haryana (INDIA)

(+91) 828 747 6144 | (+91) 875 005 0480

contact@peaar.co.in



www.aasraa.co.in