

AASRAA GLOW

FOLDING WALKING STICK WITH LED TORCH

ARL55B

The AASRAA Glow Folding Walking Stick ARL55B is a **lightweight** and **durable** walking stick designed for safety and convenience. Made from strong **aluminium**, it features a **foldable design**, **5-level height adjustment**, and a **360° pivot quad base** for **enhanced stability**. A built-in **LED torch** improves visibility in low light, while the **anti-slip base** and **ergonomic grip** ensure secure, comfortable everyday use.

★ **MRP: INR 1799.00**



ERGONOMIC AND STABLE WALKING STICK

120 KG
MAX LOAD CAPACITY

ALUMINIUM
FRAME MATERIAL

5-LEVELS
HEIGHT ADJUSTABLE

ERGONOMIC PVC
GRIP TYPE

0.4 KG
NET WEIGHT

0.5 KG
GROSS WEIGHT

15 X 9 X 90 CM
OPEN DIMENSIONS

15 X 9 X 31 CM
BOX DIMENSIONS



LED LIGHT, ERGONOMIC GRIP AND LIGHTWEIGHT

TRUCK SIZE	APPROX QTY
8'	12-16
14'	40-50
17'	70-80
22'	100-110
27'	130-140
34'	250-270
40'	290-310

LOADING CAN VARY DEPENDING ON VEHICLE INNER DIMENSIONS

SCAN FOR
PRODUCT
VIDEO &
DETAILED
SPECS



HIGHLIGHTS



SOFT LED BUTTON



STABLE SUPPORT



BUILT IN LED LIGHT



**5-LEVEL HEIGHT
ADJUSTABLE**



**PIVOT BASE WITH
360° ROTATION**



**ERGONOMIC
HANDLE**



FOLDING DESIGN



EASY PORTABLE



HIGH GRIP BASE

SPECIFICATIONS

TECHNICAL SPECIFICATIONS	
ASSEMBLY SYSTEM	Fully Assembly
FRAME MATERIAL	ALUMINIUM POWDER COATING
FRAME COLOR	BLACK POWDER COATING
FOLDABLE / FIXED	Foldable
GRIP TYPE	Ergonomic PVC Grip
ADJUSTMENT LEVEL	Height Adjustable Upto 5 Levels
FOOT TIP MATERIAL	360 Degree Pivot Anti Slip Grip

DIMENSIONS	MM	INCH
TOTAL LENGTH	150	6
TOTAL WIDTH	92	4
TOTAL HEIGHT (LOWEST TO MAX)	900-1005	35-40
GAP BETWEEN EACH HEIGHT ADJUSTMENT (RANGE)	25	1
FOLDED HEIGHT	300	12
BOX LENGTH	150	6
BOX WIDTH	90	4
BOX HEIGHT	310	20

DIMENSION	MEASUREMENT
WEIGHT OF STICK (NET WEIGHT)	0.4 KG
WEIGHT WITH PACKAGING (GROSS WEIGHT)	0.5 KG
TOTAL LEVELS OF HEIGHT ADJUSTMENT	5
LOAD CAPACITY	120 KG

SUITABLE FOR USERS



SENIORS



CHILDREN



ADULTS



WOMEN



ATHLETES



PWD



PATIENTS



PREGNANT



SHORT TERM USERS




CEREBRAL PALSRY

SUITABLE FOR ACTIVITY




FLIGHT/TRAIN



CAR



BATHROOM



OFFICE



MARKET



OUTDOOR


SUITABLE FOR SURFACE




TILE



WOOD



PAVERS




ASPHALT



CARPET



SAND



GRAVEL



GRASS



आपका भरोसा, *aasraa* का वादा

your trust, aasraa's assurance

Leading brand in India for wheelchairs (manual and electric) and other mobility and orthopaedic aids



Manufactured & Marketed By:

PEAAR HEALTHTECH LLP

Plot No 12, Sector 4
Faridabad 121004, Haryana (INDIA)

(+91) 828 747 6144 | (+91) 875 005 0480

contact@peaar.co.in



www.aasraa.co.in